Rev. DC011J011 / DC011E011

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## Safety Precautions

To help protect you and/or your property from potential damage and ensure personal safety, please read this IFU thoroughly before using the product.

### [Notational Convention]

Notation	Explanation
 CAUTION	Indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

### CAUTION

#### ■ Use only in well ventilated areas

In case of insufficient ventilation, flammable and toxic solvents can cause fire, explosion, or poisoning.

#### ■ First Aid

Skin exposure

Wash exposed area with plenty of soap and water.

Eye exposure

Open eyes as wide as possible and wash with clean water for at least 15 minutes.

Immediately call for medical attention.

Ingestion

Please wash mouth with excess water and immediately call for medical attention.

#### ■ Do not spill solvents

Spillage and leakage can cause fire, electric shock, poisoning, injury, and corrosion.

Wear appropriate protective gear when cleaning up a spill.

#### ■ Wear eye protection and protective gloves

Organic solvents and acids are harmful and should not come in direct contact with the skin.

#### ■ Handle package with care

Inappropriate handling may cause rupturing and/or spattering of the product.

#### ■ Only use this product as intended

This product is intended for *in vitro* diagnostic use for the measurement of HbA<sub>1c</sub> in blood specimens.

#### ■ Proper disposal

Dispose in accordance with local laws and regulations.

### NOTE

Keep this IFU with the product for future reference.

Европейски потребители / Evropski zákazníci / Europäische Kunden / Europäische Kunden / European customers / Ευρωπαϊκοί πελάτες / Clientes europeos / Europa klienti / Consommateurs européens / Európai ügyfelek / Clienti europei / Europas klienti / Europas klienti / Klienti europejscy / Clientes europeos / Clienti europei / Evropski kupci / Európski zákazníci / Europeiska kunder

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European Conformity

Manufacturer

Authorized  
representative in the  
European Community

Catalogue number  
/ Part number

*In vitro* diagnostic  
medical device

Consult instructions  
for use

Use by date

Batch code  
/ Lot number

Temperature limitation

For specified column lot only

Serial number  
/ Column number

Net volume  
(after reconstitution for  
lyophilized material)

## CONTENTS

1. Introduction .....	16
2. Prior to Use .....	17
3. Warnings and Precautions .....	17
4. Content .....	17
5. Part Names .....	18
6. Related Components .....	18
7. Storage .....	18
8. Specimens .....	18
9. Assay Principle .....	18
10. Installation .....	19
11. Assay Procedures .....	20
12. Precautions for Use .....	20
13. Reference Data .....	21
14. Measurement Values .....	22
15. Evaluation of Results .....	22
16. Expected Values .....	22
17. Performance Characteristics .....	22
18. Precision .....	25
19. Interferences .....	26
20. Specimen Collection and Handling .....	28

## 1. Introduction

The TSKgel G8 HSi is designed exclusively for use with the Tosoh Automated Glycohemoglobin Analyzer HLC-723G8 (referred to as HLC-723G8 in the IFU), which is based on the principle of High-Performance Liquid Chromatography assay. It is not designed for and should never be used with any other type of system.

The TSKgel G8 HSi is intended for *in vitro* diagnostic use for the measurement of hemoglobin A<sub>1c</sub> (HbA<sub>1c</sub>) in blood specimens.

Hemoglobin A<sub>1c</sub> measurements are used in the clinical management of diabetes to assess the long-term efficacy of diabetic control.

Glycohemoglobin (GHb) is a general term for complexes where the whole blood glucose is non-enzymatically bound to the  $\alpha$  or  $\beta$  chains of human hemoglobin. Within these complexes, HbA<sub>1c</sub>, which is a complex of glucose and the N terminus of the  $\beta$  chain, is the most quantitatively prevalent.

HbA<sub>1c</sub> is non-enzymatically synthesized in two steps.

Step 1: The glucose aldehyde group reacts with the free amino group on the valine comprising the N-terminus of the  $\beta$  chain to form the Schiff base (labile HbA<sub>1c</sub> (L-A<sub>1c</sub>)).

Step 2: Stable ketoamine (stable HbA<sub>1c</sub> (s-A<sub>1c</sub>)) is then formed by a reaction known as the amadori rearrangement.

Due to the fact that L-A<sub>1c</sub> (the intermediate product of this reaction) changes rapidly in response to the changes in whole blood glucose concentrations, s-A<sub>1c</sub> is now generally used to measure HbA<sub>1c</sub>. It provides the best indication of average glucose levels over the most recent 1 to 3 month period because it does not fluctuate in response to physiological factors.

In the past, it was difficult to achieve quick separation of L-A<sub>1c</sub> and s-A<sub>1c</sub> on a column, therefore sample pretreatment to remove L-A<sub>1c</sub> was necessary before

assaying s-A<sub>1c</sub>. The HLC-723G8 features quick separation of L-A<sub>1c</sub> and s-A<sub>1c</sub> on the column (1.0 min.), with no need for sample treatment when assaying s-A<sub>1c</sub>.

## 2. Prior to Use

Inspect the column package and package components for any signs of damage prior to use. If any damages are visible, contact your local Tosoh sales representative.

Confirm that the following two documents are included in the package.

- Instructions For Use 1 copy
- Glyco column inspection report 1 copy

## 3. Warnings and Precautions

- 1) This product is for *in vitro* diagnostic use only.
- 2) This product is intended for use on the Tosoh Automated Glycohemoglobin Analyzer HLC-723G8 Standard Analysis Mode only.


## 4. Content

Catalogue No.	Description	Package content
0021498	TSKgel G8 HSi	1 piece

Packed resin : Cation exchange resin  
Assay mode : Standard Analysis Mode  
Applicable analyzers : Tosoh Automated Glycohemoglobin Analyzer  
HLC-723G8

## 5. Part Names

The names of the various parts comprising the TSKgel G8 HSi are shown below.

## 6. Related Components

	Catalogue No.
Hemoglobin A1c Calibrator Set	0018767
Hemoglobin A1c Control Set	0021974
G8 Elution Buffer HSi No. 1 (S)	0021800
G8 Elution Buffer HSi No. 2 (S)	0021801
G8 Elution Buffer HSi No. 3 (S)	0021802
HSi Hemolysis & Wash Solution (L)	0018431
HSi Hemolysis & Wash Solution (LL)	0019550

## 7. Storage

Store at 4 to 15 °C in a cool dark location.

## 8. Specimens

Whole blood samples.

## 9. Assay Principle

The HLC-723G8 is based on the High-Performance Liquid Chromatography (HPLC) principle. Cation exchange columns employ the differences in ionic interactions between hemoglobin components to separate them into a total of

six fractions in a span of 1.0 min.

A step gradient elution is used to separate HbF, s-A<sub>1c</sub>, and Total A<sub>1</sub>, with three types of G8 Elution Buffer HSi (G8 Elution Buffer HSi No. 1, 2, and 3 (S)) of different salt concentrations.

## 10. Installation

Be sure to read the IFU included with the TSKgel G8 HSi, G8 Elution Buffer HSi, Hemoglobin A1c Calibrator Set, Hemoglobin A1c Control Set and HSi Hemolysis & Wash Solution as well as the HLC-723G8 Operator's Manual.

- 1) Take a column from its box and remove the end plugs from both ends. Keep the plugs in a safe place, as they will be needed when the column is stored.
- 2) Confirm that analyzer system eluent delivery has been terminated (STANDBY is displayed on the Main screen (No. 1 screen)) before opening the column oven to disconnect the tubing. Then, remove the used column.
- 3) Press the arrow key at the bottom right of the Main screen (No. 1 screen) to show the Main screen (No. 2 screen). Press the pump key on the Main screen (No. 2 screen) to start the pump and check the eluent flows through the tube disconnected from the column. Then, stop the pump by pressing the pump key again. Take special care that eluent coming from the tube does not spill onto the analyzer unit by preparing wiping paper beforehand.
- 4) Confirm the correct eluent flow direction (referring to arrow displayed on the column name plate) and connect a tube to the column inlet side. Start the pump and check that the eluent flows through the column. Then, stop the pump and connect the tube to the outlet side of the column.
- 5) Press the pump key again to start the pump. Check that the pressure level rises quickly and that there is no leakage from the tube connections.
- 6) Press the MAINT key and then the REAGENT CHANGE key. Press the COL. RESET key to reset the column counter.


## 11. Assay Procedures

Refer to the HLC-723G8 Operator's Manual for detailed instructions.

## 12. Precautions for Use


- 1) Carefully read the instructions contained in this IFU and related Instructions For Use provided with HLC-723G8, G8 Elution Buffer HSi, Hemoglobin A1c Calibrator Set, Hemoglobin A1c Control Set and HSi Hemolysis & Wash Solution.
- 2) The used column has been in contact with blood sample. Wear protective clothing (glasses, glove, mask, etc) and take sufficient care to prevent potential infection during installation and handling.
- 3) The TSKgel G8 HSi is designed exclusively for use in combination with the HLC-723G8 analyzer system, G8 Elution Buffer HSi and HSi Hemolysis & Wash Solution indicated below, never in any other combination.
  - Tosoh Automated Glycohemoglobin Analyzer HLC-723G8
  - G8 Elution Buffer HSi (No. 1, 2 and 3 (S))
  - HSi Hemolysis & Wash Solution (L ), (LL)
- 4) When replacing column, be sure to assay dummy samples to check the chromatogram results.
- 5) Always use the TSKgel G8 HSi in combination with the G8 Elution Buffer HSi of the identical lot number. The column lot number is indicated by a single uppercase alphabetical character (A, B, etc.) on the label of column box. The Elution Buffer label displays an alphabetic character corresponding to column lot number, as shown below. Note that the G8 Elution Buffer HSi No. 3 is compatible with all lots.


- 6) Care must be taken to ensure that the solutions are delivered only in the direction indicated by the arrow on the column label.
- 7) In cases in which the column is not used for more than one week, remove the column from the analyzer unit, reattach its end plugs to protect it from drying out and store in a cool dark place between 4 and 15 °C.
- 8) Please handle the column with care. Do not drop or shake the column.
- 9) If the pressure is more than 4 MPa higher than the pressure, which is indicated on the column inspection report, replace the filter. If the pressure still does not drop, please contact your local representative.
- 10) For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state and federal regulations.

### 13. Reference Data (from the HLC-723G8 Assay Examples)

#### Format 100 Control

#### Format 100 Sample

## 14. Measurement Values

Measurement values (%) indicate the percentage of each peak in relation to the Total Area (excluding the front peak (FP)). Note the minimum unit of measurement displayed is 0.1 %.

## 15. Evaluation of Results

### Quality Control

In order to monitor and evaluate the precision and analytical performance, it is recommended that control samples should be run daily. If one or more control sample value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results. Standard laboratory procedures should be followed in accordance to the strictest regulatory agency under which the laboratory operates.

## 16. Expected Values

Reference ranges (non-diabetic): HbA<sub>1c</sub> (s-A<sub>1c</sub>) 4.0-6.0 % (mean 5.0 %, SD 0.5 %)

Ref: American Diabetes Association. Standards of medical care for patients with diabetes mellitus. Diabetes Care 2000; 23 (Suppl. 1), S32-42.

The values referred to within this document have been determined with a NGSP certified method. It is known that the relationship between HbA<sub>1c</sub> results from the NGSP network (%) and the IFCC network (mmol/mol) is expressed by using the following equation (See [http://www.ifcchba1c.net/IFCC\\_08.asp](http://www.ifcchba1c.net/IFCC_08.asp)):

$$\text{NGSP} = 0.09148 \times \text{IFCC} + 2.152$$

## 17. Performance Characteristics

### 1) Linearity / Dilution

The tests were performed with Hemoglobin A1c Control diluted with various percentage factors to determine the effect on the s-A<sub>1c</sub> (%) value. Linearity was observed for the Total Area in the 500 to 2700 range, and measurements in the range of 700 to 2500 are recommended to give more reliable results.

Total Area	s-A <sub>1c</sub> (%)
324.50	9.4
652.54	9.9
939.78	9.9
1199.05	9.9
1552.53	10.0
1833.12	10.0
2104.36	10.1
2321.57	10.1
2610.74	10.1
2858.85	10.2

## 2) Recovery / Linearity


The table below shows the recovery ratios for two types of samples mixed in different percentage factors.

High s-A <sub>1c</sub> sample (ratio)	Low s-A <sub>1c</sub> sample (ratio)	Measurement value s-A <sub>1c</sub> (%)	Theoretical value s-A <sub>1c</sub> (%)	Recovery (%)
0	10	2.7	—	—
1	9	4.7	4.7	100.0
2	8	6.5	6.7	97.0
3	7	8.4	8.5	98.8
4	6	10.2	10.3	99.0
5	5	11.9	12.2	97.5
6	4	13.6	13.9	97.8
7	3	15.3	15.5	98.7
8	2	17.0	17.3	98.3
9	1	18.8	19.0	98.9
10	0	20.1	—	—


### 3) Correlation

The correlation between the s-A<sub>1c</sub> (%) values for the HLC-723G8 and HLC-723G7 analyzers and for the HLC-723G8 and HLC-723GHbV A1c2.2 analyzers are shown in the graph below. These methods are certified by the National Glycohemoglobin Standardization Program (NGSP).

The correlation between the s-A<sub>1c</sub> (%) for  
HLC-723G8 and HLC-723G7 analyzers.

The correlation between the s-A<sub>1c</sub> (%) for  
HLC-723G8 and HLC-723GHbV A1c2.2 analyzers.

## 18. Precision

### 1) Intra-assay (n=10)

The intra-assay precision (n=10) factors of the s-A<sub>ic</sub> (%) for three types of control samples are shown below.

	Low	Medium	High
Mean	5.47	7.68	9.90
SD	0.05	0.06	0.00
CV %	0.91	0.78	0.00

### 2) Inter-assay (n=10)

The inter-assay precision (n=10) factors, measured once a day for a 10-day period, of the s-A<sub>ic</sub> (%) for two types of control samples are shown below.

	Low	High
Mean	5.52	9.88
SD	0.04	0.06
CV %	0.72	0.61

## 19. Interferences


A substance is considered to interfere when the recovery of a known specimen falls outside a 10 % acceptability range.

Concentrations of interfering substances in this IFU are in units of g/L. Conversion to units of mg/dL may be made using the following equation.


$$\text{mg/dL} = \text{g/L} \times 100$$

- 1) No effect was observed below 10 g/L for glucose, 0.25 g/L for sodium cyanate and 0.20 g/L for acetaldehyde.


### Glucose

### Sodium cyanate

## Acetaldehyde

- 2) No effect was observed for the following agents up to the concentrations indicated in the table.

Free bilirubin	to 0.19 g/L
Conjugated bilirubin	to 0.20 g/L
Chyle	to 2,130 Formazine turbidity
Aspirin	to 0.50 g/L

- 3) The presence of HbD, HbS, HbC and HbE does interfere with the assay. High HbF and specific types of hemoglobin variants may interfere with the assay.
- 4) In the presence of HbE, a shoulder or small peak may be observed behind s-A<sub>1c</sub> peak.
- 5) Low concentrations of hemoglobin may cause the Total Area of the chromatogram to be < 500.
- 6) A reduced concentration of HbA<sub>1c</sub> may be observed in patients due to limited life span of red blood cells; interpret clinical results with caution.
- 7) Hemolysis does not interfere with HbA<sub>1c</sub> assay results.

## **20. Specimen Collection and Handling**

A whole blood sample is required for the assay.

No special preparation is necessary. A venous blood sample is collected aseptically. The blood samples in primary tubes, which contain EDTA, NaF, Heparin or Citric acid may be stored at 25 °C for 24 hours prior to analysis, or at 4 °C for 14 days.

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This manual may not be reprinted or copied in whole or in part without written consent of Tosoh Corporation. The contents of the manual are subject to change without notice.

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“HLC” is the registered trademark of Tosoh Corporation in Japan.

“TSKgel” is the registered trademark of Tosoh Corporation in Japan, United States, China, etc.

Printed in Japan

T1206-2